



NTA® IMPLANT

"Confidence in stability"





NTA® IMPLANT

"To improve is to change; to be perfect is to change often"
Winston Churchill

Nta Implant Company Group

Our company is engaged in providing solutions in implant dentistry by offering highest quality of dental implants. Development, manufacturing process, marketing, publicity, and trademark rights of all the products are legal and belong to our company. The products are being developed using state of the art technologies under the light of the up-to-date scientific data in the field of implantology.

We manufacture "Screw-type" and "Press-fit" implants with their accessories under the brands of Nta Implant, Nta Shorter Implant, Nta Shorter Plus Implant and Hybrid Implant. The Quality System Management Certificates belonging to us are FDA, CE, ISO 9001:2015, ISO 13485:2016 and ISO 10002:2014.

Our Solution Partners in SWITZERLAND, S.KOREA, USA and FRANCE carry out research and development, production and surface treatment for the products in their own countries. Some of our partners have more than 35 years of experience in dental market.

OUR MISSION is always to produce the most impressive products with wide product range for all dental implant indications by using effective, reliable and safe treatment methods without compromising the highest quality.

OUR VISION is to combine up-to-date developments in research and innovations in science and offer highest quality products to our customers.



Nta Implant Our Quality Certifications of Management System

ISO 9001:2015

ISO 13485:2016

ISO 10002:2014



6
PAGE

15
PAGE

18
PAGE

32
PAGE

INFORMATION

IMPLANT


ABUTMENT

SURGICAL KIT

DIAMETER : Ø 3.4 / Ø 3.7 / Ø 4.2 / Ø 4.6 / Ø 5.1 / Ø 6.0
 LENGTH : 6.6 / 7.5 / 8.5 / 10 / 11.5 / 13
 SURFACE : SLA

Platform Design

Smaller platform compared to many designs allows platform switching.



Upper Thread

Taper thread on the straight body provides even distribution of the load and exerts the stress on the implant. Lower depth of upper threads helps to promote vascularisation on the cortical bone area and therefore reduces the time required for healing and new bone formation.

Apex

The rounded apex may reduce the risk of tissue damage, especially when placing the implant into maxillary sinus sites.

Conical Connection With Hex 11° Index


Conical connection provides bacterial seal and prevents microleakage issues at implant-abutment interface. Internal hex index ensures proper and easy seating of the abutment.

Confidence in Stability



Single Pitch Macro Thread

- Synchronized thread
- Preferred design for SLA surface
- High fixture strength



Open Thread

- Prevents bone necrosis
- Reduces stress concentration

Corkscrew Thread

- Superior self-threading
- Ensures implant path
- Easy path change
- High primary stability
- Even stress distribution
- Decrease sensitivity on drill size

Helix Cutting Edge

- Superior self-threading
- Easy path change

Apical Design

- Self-drilling ability
- Increased initial stability in the presence of insufficient bone quality and volume

“Confidence in Stability!”




The Nta Sla Implant Surface

The SLA surface is produced by a large grit sand-blasting process with corundum particles that leads to a macro roughness on the titanium surface. This is followed by a strong acid-etching bath with a HCl at elevated temperature for several minutes. This produces the fine 2–4 µm micro pits superimposed on the rough-blasted surface.

The surface is not micro porous and therefore provides no enclosed volumes to reduce vulnerability to bacteria.

The chemical composition of the SLA structure was found to be titanium oxide (TiO_2) using X-ray photoelectron spectroscopy.

This method analyses the first few atomic layers of the surface, and thus the chemical composition of the material which is in direct contact and interacts with tissue fluids and cells.



No-Mount System

Surgery time has been shortened by removing the fixture mount from the package and by using a fixture driver for which the mount is not necessary. However, fixture mount can be provided separately according to the preference of the clinicians.

Surgical Drill

Minimize the use of drill extension by making the full length of drill to be 33mm. It can be used conveniently with its 1mm counterbore function in cases where the bone is not flat or when the placement of the implant below the bone level is required.



Conical Connection With Internal Index

The implant has a conical connection, which provides bacterial seal, eliminates the micromovement and micoleakage issues. The bone structure is therefore preserved and abutment screw loosening is prevented. Internal hex index on the other hand, ensures proper and precise seating of the abutment.



Platform Switching


Implant-abutment junction is placed closer to the center of the implant and away from implant-bone interface. Together with advantages of conical connection, crestal bone loss is a rare phenomenon in NTA implants.

Flange Cutting

Cutting properties of the threads minimize the trauma caused during implant insertion.


Cutting Edge

When using the drills that contain three sharp edges, the preparation of the implant osteotomy is easy, the frictional heat formation is low, post-operative pain is limited and the drilling time is quite short.



Root Form

Root form of the implant ensures easier placement, high initial stability, superior stress distribution, less risk of damaging neighboring roots and low risk for accidental perforation of concave surfaces.



Non Fixture

If the neck area is exposed, easy cleaning of the area it provided so that the possibility of having perimucositis and periimplantitis is reduced.

NTA COMPLETE GUIDE



NTA PACKING SYSTEM

New Ampoule User's Guide



Open the package.

Remove the back of blister (sterility is not guaranteed if the seal is broken).

Keep the small cover upward.


Remove the cover by gently rotating and pulling.

Insert the driver into the fixture and make sure that the driver is fully seated.


Closing Screw
Gold metallic color
Ø 3.4 narrowClosing Screw
Silver metallic color
Ø 3.7 - 4.2 - 4.6 - 5.1 - 6.0 regular

NTA SURGICAL KIT PROTOKOL


BONE	MEDIUM & HARD (D1 + D2+ D3)	SOFT (D4)
F	3.4	2.0 + 3.4
I	3.7	2.0 + 3.4 + 3.7
X	4.2	2.0 + 3.4 + 3.7 + 4.2
T	4.6	2.0 + 3.4 + 3.7 + 4.2 + 4.6
U	5.1	2.0 + 3.4 + 3.7 + 4.2 + 4.6 + 5.1
R	6.0	2.0 + 3.4 + 3.7 + 4.2 + 4.6 + 5.1 + 6.0
E		2.0 + 3.4 + 3.7 + 4.2 + 4.6 + 5.1 + 6.0




NTA PROSTHETIC FLOW DIAGRAMS FOR OCTA SYSTEM




NTA PROSTHETIC FLOW DIAGRAMS FOR MULTI (ALL ON FOUR) SYSTEM



NTA PROSTHETIC FLOW DIAGRAMS FOR NTA SYSTEM
OVERDENTURE RETAINED RESTORATION : N-LOC




Internal Sub-merged System



"CONFIDENCE IN STABILITY"

Narrow

$\varnothing 3.4$
Morse Tapered Connection 11°



Straight Body

As implant length gets longer, the straight body gets longer. For this reason, long implants can be placed steadily without causing unexpected changes at the time of the surgery same as straight type of implants.

Morse Tapered Connection 11°

Implant and abutment excellent connection is guaranteed. 11° degree morse tapered connection is reducing stress on the fixture and decrease bone loss.

Open Thread



In the course of implant placement, open thread gives opportunity to place fixture deep without necessity of additional drilling by minimising resistance.

Cutting Edge

The sharp cutting edge maximized self-tapping.

Regular

$\varnothing 3.7 - \varnothing 4.2 - \varnothing 4.6 - \varnothing 5.1 - \varnothing 6.0$
Morse Tapered Connection 11°



Fixture

$\varnothing 3.4 - \varnothing 3.7$

TF TYPE $\varnothing 3.4 \rightarrow$ NARROW 2.1	Length	Product Code	Screw
	8.5	NTF 3408 S	
	10	NTF 3410 S	
	11.5	NTF 3411 S	M 1.6 Abutment Screw M 1.6
	13	NTF 3413 S	

TF TYPE $\varnothing 3.7 \rightarrow$ REGULAR 2.5	Length	Product Code	Screw
	7.5	NTF 3707 S	
	8.5	NTF 3708 S	
	10	NTF 3710 S	M 2.0 Abutment Screw M 2.0
	11.5	NTF 3711 S	
	13	NTF 3713 S	

NM TYPE $\varnothing 4.2 \rightarrow$ REGULAR 2.5	Length	Product Code	Screw
	7.5	FNM 4007 S	
	8.5	FNM 4008 S	
	10	FNM 4010 S	M 2.0 Abutment Screw M 2.0
	11.5	FNM 4011 S	
	13	FNM 4013 S	

NM TYPE $\varnothing 4.6 \rightarrow$ REGULAR 2.5	Length	Product Code	Screw
	6.6	FNM 4506 S	
	7.5	FNM 4507 S	
	8.5	FNM 4508 S	
	10	FNM 4510 S	M 2.0 Abutment Screw M 2.0
	11.5	FNM 4511 S	
	13	FNM 4513 S	

NM TYPE $\varnothing 5.1 \rightarrow$ REGULAR 2.5	Length	Product Code	Screw
	6.6	FNM 5006 S	
	7.5	FNM 5007 S	
	8.5	FNM 5008 S	
	10	FNM 5010 S	M 2.0 Abutment Screw M 2.0
	11.5	FNM 5011 S	
	13	FNM 5013 S	

NM TYPE $\varnothing 6.0 \rightarrow$ REGULAR 2.5	Length	Product Code	Screw
	6.6	FNM 6006 S	
	7.5	FNM 6007 S	
	8.5	FNM 6008 S	M 2.0 Abutment Screw M 2.0


No Mount System: Uses a special fixture driver, which is more durable and has easy separation and connection properties.

Surface: Sha
Recommended Placement Torque: Less than 40 Ncm.
Tapered angle for safer maxillary placement.





**Esthetic!
Abutment**

*Esthetic Abutments conical connection;
Prefabricated titanium abutments available in straight and angled;
and varying collar heights.*



NTA HEALING ABUTMENT




- ✓ Hex Driver: 1.2 mm.
- ✓ Fastening Torque: 5-8 Ncm.

Diameter	*G/H	Height	Product Code
Narrow Ø 4.15	2	1	HA 403 M
	3	1.5	HA 405 M
	4	2.5	HA 407 M
	5	2.5	HA 409 M
Narrow Ø 4.5	2	1	HA 453 M
	3	1.5	HA 455 M
	4	2.5	HA 457 M
	5	2.5	HA 459 M
Narrow Ø 5.5	2	1	HA 553 M
	3	1.5	HA 555 M
	4	2.5	HA 557 M
	5	2.5	HA 559 M
Regular Ø 4.15	1	1	HA 402 S
	2	1	HA 403 S
	3	1.5	HA 405 S
	4	2.5	HA 407 S
	5	2.5	HA 409 S
Regular Ø 4.5	1	1	HA 452 S
	2	1	HA 453 S
	3	1.5	HA 455 S
	4	2.5	HA 457 S
	5	2.5	HA 459 S
Regular Ø 5.5	1	1	HA 552 S
	2	1	HA 553 S
	3	1.5	HA 555 S
	4	2.5	HA 557 S
	5	2.5	HA 559 S
Regular Ø 6.5	1	1	HA 652 S
	2	1	HA 653 S
	3	1.5	HA 655 S
	4	2.5	HA 657 S

*GH : Gingival Height




+ NTA TEMPORARY ABUTMENT



Fixture	Diameter	*G/H	Height	Product Code
Narrow	Ø 4.5	1	8	NTTAN 4510 H
		2	8	NTTAN 4520 H
		3	8	NTTAN 4530 H
Regular	Ø 4.5	1	8	NTTAR 4510 H
		2	8	NTTAR 4520 H
		3	8	NTTAR 4530 H

✓ Short screw should be used with short implants.

+ NTA TEMPORARY ABUTMENT - TITANIUM



Fixture	Diameter	*G/H	Height	Product Code
Narrow	Ø 4.5	1	10	NTTTAN 4510 H
Regular	Ø 4.5	1	10	NTTTAR 4510 H



✓ Short screw should be used with short implants.

+ NTA SOLID ABUTMENT



Fixture	Diameter	*G/H	Height	Product Code
Narrow	Ø 4.5	1	4.0	NTSA 30514 M
		2		NTSA 30524 M
		3		NTSA 30534 M
		4		NTSA 30544 M
		1	5.5	NTSA 30515 M
		2		NTSA 30525 M
		3		NTSA 30535 M
		4		NTSA 30545 M
		1	7.0	NTSA 30517 M
		2		NTSA 30527 M
		3		NTSA 30537 M
		4		NTSA 30547 M
		1	7.0	NTSA 30617 M
		2		NTSA 30627 M
		3		NTSA 30637 M
		4		NTSA 30647 M
Ø 5.5	Ø 5.5	1	4.0	NTSA 65714
		2		NTSA 65724
		3		NTSA 65734
		4		NTSA 65744
		1	5.5	NTSA 65715
		2		NTSA 65725
		3		NTSA 65735
		4		NTSA 65745
		5	5.5	NTSA 65755
		6		NTSA 65765
		1	7.0	NTSA 65717
		2		NTSA 65727
		3		NTSA 65737

+ NTA SOLID ABUTMENT



Fixture	Diameter	*G/H	Height	Product Code
Regular	Ø 4.5	1	4.0	NTSA 45514
		2		NTSA 45524
		3		NTSA 45534
		4		NTSA 45544
		1	5.5	NTSA 45515
		2		NTSA 45525
		3		NTSA 45535
		4		NTSA 45545
		5	7.0	NTSA 45555
		6		NTSA 45565
		1		NTSA 45517
		2		NTSA 45527
		3		NTSA 45537
		4		NTSA 45547
		1	4.0	NTSA 55614
Ø 5.5	Ø 5.5	2		NTSA 55624
		3		NTSA 55634
		4		NTSA 55644
		1		NTSA 55615
		2		NTSA 55625
		3		NTSA 55635
		4		NTSA 55645
		5		NTSA 55655
		6		NTSA 55665
		1	7.0	NTSA 55617
		2		NTSA 55627
		3		NTSA 55637
		4		NTSA 55647
Ø 6.5	Ø 6.5	1	4.0	NTSA 65714
		2		NTSA 65724
		3		NTSA 65734
		4		NTSA 65744
		1	5.5	NTSA 65715
		2		NTSA 65725
		3		NTSA 65735
		4		NTSA 65745
		5	5.5	NTSA 65755
		6		NTSA 65765
		1	7.0	NTSA 65717
		2		NTSA 65727
		3		NTSA 65737

- ✓ An abutment and a set screw, used for fixed prostheses.
- ✓ Hex driver: 1.2mm.
- ✓ Fastening Torque: 20Nm (Narrow) - 30Nm (Regular)
- ✓ Protect cap: Optional.

*G/H : Gingival Height

NTA CEMENTED ABUTMENT - HEX


Fixture	Diameter	*G/H	Height	Product Code	Screw
Narrow	Ø 4.15	1	4.0	NTCA 20514 H	M 1.6 AAHS 1610
		2		NTCA 20524 H	
		3		NTCA 20534 H	
		4		NTCA 20544 H	
		1	5.5	NTCA 20515 H	
		2		NTCA 20525 H	
		3		NTCA 20535 H	
		4		NTCA 20545 H	
		5	7.0	NTCA 20555 H	
		6		NTCA 20565 H	
		1		NTCA 20517 H	
		2		NTCA 20527 H	
		3		NTCA 20537 H	
		4		NTCA 20547 H	
Regular	Ø 4.5	1	4.0	NTCA 30514 H	M 1.6 AAHS 1610
		2		NTCA 30524 H	
		3		NTCA 30534 H	
		4		NTCA 30544 H	
		1	5.5	NTCA 30515 H	
		2		NTCA 30525 H	
		3		NTCA 30535 H	
		4		NTCA 30545 H	
		5	7.0	NTCA 30555 H	
		6		NTCA 30565 H	
		1		NTCA 30517 H	
		2		NTCA 30527 H	
		3		NTCA 30537 H	
		4		NTCA 30547 H	
Regular	Ø 4.15	1	4.0	NTCA 40514 H	M 2.0 AAHS 2010
		2		NTCA 40524 H	
		3		NTCA 40534 H	
		4		NTCA 40544 H	
		1	5.5	NTCA 40515 H	
		2		NTCA 40525 H	
		3		NTCA 40535 H	
		4		NTCA 40545 H	
		5	7.0	NTCA 40555 H	
		6		NTCA 40565 H	
		1		NTCA 40517 H	
		2		NTCA 40527 H	
		3		NTCA 40537 H	
		4		NTCA 40547 H	

✓ Abutment for both single unit prosthesis and multiunit fixed partial restorations.
 ✓ Hex driver: 1.2mm.
 ✓ Fastening Torque: 20Ncm (Narrow) - 30Ncm (Regular)
 ✓ Protect cap: Optional.




NTA CEMENTED ABUTMENT - HEX


Fixture	Diameter	*G/H	Height	Product Code	Screw
Regular	Ø 4.5	1	4.0	NTCA 45514 H	M 2.0 AAHS 2010
		2		NTCA 45524 H	
		3		NTCA 45534 H	
		4		NTCA 45544 H	
		1	5.5	NTCA 45515 H	
		2		NTCA 45525 H	
		3		NTCA 45535 H	
		4		NTCA 45545 H	
		5	7.0	NTCA 45555 H	
		6		NTCA 45565 H	
		1		NTCA 45517 H	
		2		NTCA 45527 H	
		3		NTCA 45537 H	
		4		NTCA 45547 H	
Regular	Ø 5.5	1	4.0	NTCA 55614 H	M 2.0 AAHS 2010
		2		NTCA 55624 H	
		3		NTCA 55634 H	
		4		NTCA 55644 H	
		1	5.5	NTCA 55615 H	
		2		NTCA 55625 H	
		3		NTCA 55635 H	
		4		NTCA 55645 H	
		5	7.0	NTCA 55655 H	
		6		NTCA 55665 H	
		1		NTCA 55617 H	
		2		NTCA 55627 H	
		3		NTCA 55637 H	
		4		NTCA 55647 H	
Regular	Ø 6.5	1	4.0	NTCA 65714 H	M 2.0 AAHS 2010
		2		NTCA 65724 H	
		3		NTCA 65734 H	
		4		NTCA 65744 H	
		1	5.5	NTCA 65715 H	
		2		NTCA 65725 H	
		3		NTCA 65735 H	
		4		NTCA 65745 H	
		5	7.0	NTCA 65755 H	
		6		NTCA 65765 H	
		1		NTCA 65717 H	
		2		NTCA 65727 H	
		3		NTCA 65737 H	
		4		NTCA 65747 H	

✓ Abutment for both single unit prosthesis and multiunit fixed partial restorations.
 ✓ Hex driver: 1.2mm.
 ✓ Fastening Torque: 20Ncm (Narrow) - 30Ncm (Regular)
 ✓ Protect cap: Optional.

*GH : Gingival Height


+ NTA CEMENTED ABUTMENT - NON-HEX



Fixture	Diameter	*G/H	Height	Product Code	Screw
Regular	Ø 4.5	1	4.0	NTCA 45514 NH	M 2.0 AAHS 2010
		2		NTCA 45524 NH	
		3		NTCA 45534 NH	
		4		NTCA 45544 NH	
	Ø 4.5	1	5.5	NTCA 45515 NH	
		2		NTCA 45525 NH	
		3		NTCA 45535 NH	
		4		NTCA 45545 NH	
	Ø 5.5	5	4.0	NTCA 45555 NH	
		6		NTCA 45565 NH	
		1	7.0	NTCA 45517 NH	
		2		NTCA 45527 NH	
		3		NTCA 45537 NH	
		4		NTCA 45547 NH	
	Ø 6.5	1	4.0	NTCA 55614 NH	
		2		NTCA 55624 NH	
		3		NTCA 55634 NH	
		4		NTCA 55644 NH	
		1	5.5	NTCA 55615 NH	
		2		NTCA 55625 NH	
		3		NTCA 55635 NH	
		4		NTCA 55645 NH	
		5	7.0	NTCA 55655 NH	
		6		NTCA 55665 NH	
		1		NTCA 55617 NH	
		2		NTCA 55627 NH	
		3		NTCA 55637 NH	
		4		NTCA 55647 NH	
		1	4.0	NTCA 65714 NH	
		2		NTCA 65724 NH	
		3		NTCA 65734 NH	
		4		NTCA 65744 NH	
		1	5.5	NTCA 65715 NH	
		2		NTCA 65725 NH	
		3		NTCA 65735 NH	
		4		NTCA 65745 NH	
		5	7.0	NTCA 65755 NH	
		6		NTCA 65765 NH	
		1		NTCA 65717 NH	
		2		NTCA 65727 NH	
		3		NTCA 65737 NH	
		4		NTCA 65747 NH	





✓ Abutment for both single unit prosthesis and multiunit fixed partial restorations.
✓ Hex driver: 1.2mm.
✓ Fastening Torque: 20Ncm (Narrow) - 30Ncm (Regular)
✓ Protect cap: Optional.

+ NTA MILLING ABUTMENT



Fixture	Diameter	*G/H	Height	Product Code-Hex	Product Code Non-Hex
Regular	Ø 4.5	2	12	NTMA 452 MH	NTMA 452 MNH
	Ø 4.5	2		NTMA 452 H	NTMA 452 NH
	Ø 5.5	2		NTMA 552 H	NTMA 552 NH
	Ø 5.5	3		NTMA 553 H	NTMA 553 NH
	Ø 6.5	2		NTMA 652 H	NTMA 652 NH
	Ø 6.5	4		NTMA 654 H	NTMA 654 NH

+ NTA ANGLED ABUTMENT-HEX




Fixture	Diameter	*G/H	Height	Angle	Product Code	Screw
Narrow	Ø 4.5	1	7.0	7°	NTAA 45107 MH	M 1.6 AAHS 1608A
		2		7°	NTAA 45207 MH	
		4		7°	NTAA 45407 MH	
		1		15°	NTAA 45115 MH	
		2		15°	NTAA 45215 MH	
		4		15°	NTAA 45415 MH	
		1		25°	NTAA 45125 MH	
		2		25°	NTAA 45225 MH	
Regular	Ø 5.5	4		25°	NTAA 45425 MH	
		2	7.0	7°	NTAA 55207MH	
		4		7°	NTAA 55407MH	
		2		15°	NTAA 55215 MH	
		4		15°	NTAA 55415 MH	
		2		25°	NTAA 55225 MH	
		4		25°	NTAA 55425 MH	
		1	7.0	7°	NTAA 50107 SH	M 2.0 AAHS 2008B
Ø 5.0	Ø 5.0	2		7°	NTAA 50207 SH	
		4		7°	NTAA 50407 SH	
		1		15°	NTAA 50115 SH	
		2		15°	NTAA 50215 SH	
		4		15°	NTAA 50415 SH	
		1		25°	NTAA 50125 SH	
		2		25°	NTAA 50225 SH	
		4		25°	NTAA 50425 SH	
Ø 6.0	Ø 6.0	2	7.0	7°	NTAA 60207 SH	
		4		7°	NTAA 60407 SH	
		2		15°	NTAA 60215 SH	
		4		15°	NTAA 60415 SH	
		2		25°	NTAA 60225 SH	
		4		25°	NTAA 60425 SH	

✓ Fastening Torque: 20Ncm (Narrow) - 30Ncm (Regular)

*GH : Gingival Height


+ NTA ANGLED ABUTMENT - NON-HEX




Fixture	Diameter	*G/H	Height	Angle	Product Code	Screw
Regular	Ø 5.0	2	7.0	15°	NTAA 50215 NH	M 2.0 SCA 2008 B
		4		15°	NTAA 50415 NH	
		2		25°	NTAA 50225 NH	
		4		25°	NTAA 50425 NH	
	Ø 6.0	2	7.0	15°	NTAA 60215 NH	
		4		15°	NTAA 60415 NH	
		2		25°	NTAA 60225 NH	
		4		25°	NTAA 60425 NH	

✓ Fastening Torque: 30Ncm (Regular)

+ NTA ESTHETIC ANGLED ABUTMENT



Fixture	Diameter	*G/H	Angle	Height	Product Code
Narrow	4.0	1.5	7°	7.0	NTEA 274015
		2	7°		NTEA 274020
		3	7°		NTEA 274030
		4	7°		NTEA 274040
		1.5	17°		NTEA 24015
	4.5	2	17°		NTEA 24020
		3	17°		NTEA 24030
		4	17°		NTEA 24040
		1.5	7°		NTEA 274515
		2	7°		NTEA 274520



Fixture	Diameter	*G/H	Angle	Height	Product Code
Regular	4.5	1.5	7°	7.0	NTEA 574515
		2	7°		NTEA 574520
		3	7°		NTEA 574530
		4	7°		NTEA 574540
		1.5	17°		NTEA 54515
	4.5	2	17°		NTEA 54520
		3	17°		NTEA 54530
		4	17°		NTEA 54540
		1.5	7°		NTEA 575515
		2	7°		NTEA 575520


✓ Fastening Torque: 20Ncm (Narrow) - 30Ncm (Regular)

+ NTA UCLA ABUTMENT




Fixture	Product Code	Screw
Narrow	UCA2107 M	AAHS 1610
Regular	UCA2507	SURO-S

Short screw should be used with short implants.



Esthetic Abutments conical connection;
Prefabricated titanium abutments available in straight and angled ; and varying collar heights.

✓ Fastening Torque: 20Ncm (Narrow) - 30Ncm (Regular)



NTA OCTA ABUTMENT



Fixture	Diameter	*G/H	Product Code	Low Plastic Cylinder
Narrow	$\varnothing 4.0$	1	OAA 481 M	POI4810N (Non-Hex) POI4810H (Hex)
		2	OAA 482 M	
		3	OAA 483 M	
		4	OAA 484 M	
Regular	$\varnothing 4.0$	1	OAA 481	
		2	OAA 482	
		3	OAA 483	
		4	OAA 484	

✓ Fastening Torque: 20Ncm (Narrow) - 30Ncm (Regular)

NTA OCTA ABUTMENT PROTECT CAP

Diameter	Height	Product Code
$\varnothing 5.0$	4	OAPC 4840

✓ Torx Driver: 1.7 mm.
✓ Fastening Torque: 20Ncm.
✓ Packing Contents: Healing cap + Cylinder screw.

NTA OCTA ABUTMENT PICK-UP IMPRESSION COPING

Diameter	Height	Product Code
$\varnothing 5.5$	10	OAIC 4810

✓ Hex driver: 1.2 mm.
✓ Minimizes diskortion.
✓ Packing Contents: impression coping + Guided pin.

NTA OCTA ABUTMENT LAB FIXTURE ANALOG

Diameter	Height	Product Code
$\varnothing 5.0$	14	OALFA 4814




NTA MULTI ABUTMENT ALL ON FOUR SYSTEM



Fixture	Diameter	*G/H	Product Code
Narrow	$\varnothing 4.8$	1	NTMA(N) 5010
		2	NTMA(N) 5020
		3	NTMA(N) 5030
		4	NTMA(N) 5040
Regular	$\varnothing 4.8$	1	NTMA(R) 5010
		2	NTMA(R) 5020
		3	NTMA(R) 5030
		4	NTMA(R) 5040

✓ Fastening Torque: 20Ncm (Narrow) - 30Ncm (Regular)

NTA MULTI ANGLED ABUTMENT ALL ON FOUR SYSTEM



Fixture	Diameter	*G/H	Angle	Product Code	Screw
Narrow	$\varnothing 4.8$	2	17°	NT17MAN 4820	BMANS200
		3	17°	NT17MAN 4830	
		4	17°	NT17MAN 4840	
		3	30°	NT30MAN 4830	
		4	30°	NT30MAN 4840	
		5	30°	NT30MAN 4850	
Regular	$\varnothing 4.8$	2	17°	NT17MAR 4820	BMARS200
		3	17°	NT17MAR 4830	
		4	17°	NT17MAR 4840	
		3	30°	NT30MAR 4830	
		4	30°	NT30MAR 4840	
		5	30°	NT30MAR 4850	

✓ Fastening Torque: 20Ncm (Narrow) - 30Ncm (Regular)

NTA MULTI ABUTMENT LOW TRANSFER IMPRESSION COPING



Diameter	Height	Product Code
4.8	8.0	METIC 4880

NTA MULTI ABUTMENT LOW PLASTIC CYLINDER

Product Code	Diameter
MEPC 4846	4.8

NTA MULTI ABUTMENT OUTER DRIVER

Product Code
MAOD



NTA MULTI ABUTMENT LOW TEMPORARY CYLINDER

Product Code
METC 48

NTA MULTI ABUTMENT LOW HEALING CAP


Diameter	Height	Product Code
4.8	4.6	MEHC 4880



Product Code
MELA 300

Product Code
MEOTIC 4811


NTA BALL ABUTMENT



Fixture	Diameter	*G/H	Product Code
Narrow	$\varnothing 4.5$	1	NTBA 4510 M
		2	NTBA 4520 M
		3	NTBA 4530 M
		4	NTBA 4540 M
		5	NTBA 4550 M
		6	NTBA 4560 M
Regular	$\varnothing 4.5$	1	NTBA 4510
		2	NTBA 4520
		3	NTBA 4530
		4	NTBA 4540
		5	NTBA 4550
		6	NTBA 4560

✓ Abutments for removable over-dentures.

✓ Fastening Torque: 20Ncm (Narrow) - 30Ncm (Regular)




NTA T-BASE ABUTMENT

Fixture	Diameter	*G/H	Product Code
Narrow	Ø 4.5	1	NTBA NR
Regular	Ø 4.5	1	NTBA RG

✓ Fastening Torque: 20Ncm (Narrow) - 30Ncm (Regular)

NTA LOC



Fixture	*G/H	Angel	Code
Narrow	1	0°	NT 401 N
	1,5	15°	NT 41515 N
	2	0°	NT 402 N
	3	0°	NT 403 N
	3	15°	NT 40315 N
	4	0°	NT 404 N
	5	0°	NT 405 N
	6	0°	NT 406 N


✓ Pick-up type (to be used with a custom tray).



Fixture	*G/H	Angel	Code
Regular	1	0°	NT 401 R
	1,5	15°	NT 41515 R
	2	0°	NT 402 R
	3	0°	NT 403 R
	3	15°	NT 40315 R
	4	0°	NT 404 R
	5	0°	NT 405 R
	6	0°	NT 406 R




NTA TRANSFER IMPRESSION COPING



Fixture	Diameter	Height	Product Code	Pin
Narrow	Ø 4.0	11 (Short)	TIC 4011 M	TICP 1617
		14 (Long)	TIC 4014 M	TICP 1620
	Ø 5.0	11 (Short)	TIC 5011 M	TICP 1617
		14 (Long)	TIC 5014 M	TICP 1620
Regular	Ø 4.0	11 (Short)	TIC 4011 S	TICP 2017
		14 (Long)	TIC 4014 S	TICP 2020
	Ø 5.0	11 (Short)	TIC 5011 S	TICP 2017
		14 (Long)	TIC 5014 S	TICP 2020
	Ø 6.0	11 (Short)	TIC 6011 S	TICP 2017
		14 (Long)	TIC 6014 S	TICP 2020

✓ Impression coping for close tray.
✓ Packing Unit: Impression coping + guide pin.

✓ Hex Type.
✓ Available in short and long sizes.




NTA PICK-UP IMPRESSION COPING - HEX



Fixture	Diameter	Height	Product Code	Pin
Narrow	Ø 4.5	10 (Short)	PIC 3410 MH	PICP 1613
		13 (Long)	PIC 3413 MH	PICP 1617
Regular	Ø 4.5	10 (Short)	PIC 4810 SH	PICP 2013
		13 (Long)	PIC 4813 SH	PICP 2017
	Ø 5.5	10 (Short)	PIC 5810 SH	PICP 2013

✓ Impression coping for custom tray (to be used with hex abutments).
✓ Available in short and long sizes.

NTA PICK-UP IMPRESSION COPING - NON-HEX



Fixture	Diameter	Height	Product Code	Pin
Narrow	Ø 4.5	10 (Short)	PIC 3410 MNH	PICP 1613
		13 (Long)	PIC 3413 MNH	PICP 1617
Regular	Ø 4.5	10 (Short)	PIC 4810 SNH	PICP 2013
		13 (Long)	PIC 4813 SNH	PICP 2017
	Ø 5.5	10 (Short)	PIC 5810 SNH	PICP 2013

✓ Impression coping for custom tray (to be used with non-hex abutments).
✓ Available in short and long sizes.

NTA LAB FIXTURE ANALOG



Fixture	Diameter	Product Code
Narrow	Ø 3.4	LFA 3412 M
Regular	Ø 3.7, Ø 4.2, Ø 4.6, Ø 5.1 Ø 6.0	LFA 4012

✓ Fixture analogs for the lab (to be used with appropriate abutments)

NTA SOLID LAB ANALOG



Diameter	Height	Product Code
Ø 4.5	4	SLA 4504
	5.5	SLA 4505
Ø 5.5	5.5	SLA 5505
	7	SLA 5507
Ø 6.5	4	SLA 6504
	5.5	SLA 6505
Ø 6.5	7	SLA 6507






✓ Abutment analog (to be used with abutment level impression technique).

NTA IMPRESSION PLASTIC COPING




Diameter	Product Code
Ø 4.15	IPC 4515
Ø 4.5	IPC 4510
Ø 5.5	IPC 5510
Ø 6.5	IPC 6510

✓ To be used both with cemented and solid abutments.

NTA SURGICAL KIT PROTOKOL




Surgical Tool Maintenance Process



Multi Abutment!
All On Four System

*Ease of use
Intelligent design
Restorative flexibility*



Multi Abutment



EC REP

AZ & PARTNER AG

Dorfchärrn CH-6243 Egolzwil, SWITZERLAND

Tel: +41 41 982 02 02 Fax: +41 41 982 03 03

www.ntaimplant.com - info@ntaimplant.com

NTA Implant Ltd.

İsmetpaşa Cad. 43-3 Antalya / TURKEY

Tel: +90 242 248 19 96 Fax: +90 242 244 02 14